EN: User Manual

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SV: Bruksanvisning

DA: Brugervejledning

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Head office:
 GRANBERG AS
 Bjoavegen 1442, NO-5584 Bjoa, Norway
 Tel +47 53 77 53 00, e-mail post@granberg.no

Swedish office / European warehouse:
 GRANBERG SVERIGE AB
 Schubergvägen 20, 311 74 Falkenberg, Sweden
 Tel +46 (0)346 124 25, e-mail post@granberg-ab.se

www.granberggloves.com

User Manual issue date: 27.09.2023

Explanation of symbols and pictograms used

- Consult instructions for use
- Corrugated cardboard
- Manufacturer
- Non-corrugated cardboard
- Date of manufacture
- Paper
- Expiry date
- Low density polyethylene
- Lot number
- Medical Device
- Fragile, handle with care
- Unique Device Identifier
- Article number
- Keep away from sunlight
- Keep dry
- Temperature limit
- Non-sterile
- Do not reuse
- Suitable for contact with foodstuffs. Check the Food Declaration of Compliance.

ENGLISH

210.0022E – GRANBERG SURGICAL DISPOSABLE FACE MASKS TYPE IIR

Description

Disposable surgical face masks Type IIR are made of non-woven polypropylene. Material in 3-layers with elastic loop for ears. The face masks are 100% latex free. >99% bacterial filtration efficiency.

The face masks are made of following material:

- Outer layer:* Polypropylene.
- Second layer:* Meltblown fabric.
- Inner material:* Polypropylene.
- Ear loops:* Polyurethane.
- Nose clip material:* Polypropylene and steel.

Medical Device Regulation

This product is classified as a Class I medical device and meets the provisions of the Medical Device Regulation (EU) 2017/745. Non-sterile, single use.

Intended purpose: medical face masks are designed to prevent transmission of infective agents, blood and/or bodily fluids from the person wearing them (including but not limited to healthcare personnel and patients) onto others.

Performance

This product meets the requirements of the European Standard EN 14683:2019+AC:2019: Surgical face masks, type IIR.

Test	Type I *	Type II	Type IIR
Bacterial filtration efficiency (BFE), (%)	≥ 95	≥ 98	≥ 98
Differential pressure (Pa/cm ²)	< 40	< 40	< 60
Splash resistance pressure (kPa)	Not required	Not required	≥ 16,0
Microbial cleanliness (cfu/g)	≤ 30	≤ 30	≤ 30

* Type I medical face masks should only be used for patients and other persons to reduce the risk of spread of infections particularly in epidemic or pandemic situations. Type I masks are not intended for use by healthcare professionals in an operating room or in other medical settings with similar requirements.

Face masks in compliance with the Regulation (EC) 1935/2004 and according to Commission Regulation (EU) No. 10/2011 suitable for handling all types of foods for a short period of time.

Material Performance

The performance of the face masks quoted above is based on laboratory data and may not reflect the actual duration of protection in the workplace due to other factors influencing the performance such as temperature, mechanical wear etc.

User guide

1. Wash your hands thoroughly before you put on the mask.
2. Hold the mask by its edges in front of you, with the white side inwards and the nose clip upwards. Pull a rubber band behind each ear.
3. Pinch and stroke your fingers along the nose clip a couple of times to form it according to the shape of your nose ridge, preventing air from entering between the mask and your nose.
4. Pinch the nose clip while pulling the lower part of the face mask below your chin to achieve a better and tighter fit.
5. Do not touch the mask while wearing it, and only hold it by the rubber bands when removing it. Discard the mask in a proper manner.

Do not use the mask for more than 4 hours. Remove it if breathing becomes difficult, it becomes damp or the fit of the mask becomes loose.

Storage

This product should be stored in a dry, ventilated environment where the temperature is between +5°C to +30°C. Keep away from direct sunlight, fire and flammable substances. Preferably stored in original packaging.

When stored as recommended, the chemical properties will remain unchanged for up to 2 years from the date of manufacture. Lifetime cannot be specified as it depends on its use. It is the responsibility of the user to determine whether the face mask is suitable for intended use.

Reporting of serious incidents

In case of any serious incident has occurred in relation to the device, please report to the manufacturer and the competent Authority.

Use limitations and Disposal

- The used raw materials are chosen to minimize the risk of allergy and hypersensitivity, but a guarantee for allergic reactions cannot be given. If any allergic reaction occurs, seek medical advice or visit a doctor immediately.
- EU Declaration of Conformity is accessible at www.granberg.no/catalog/210.0022E
- For single use only. Do not re-use.
- Before usage, inspect face mask for any defect or imperfections.
- Dispose according to hospital policy and local regulation.
- The manufacturer cannot be held responsible for improper use of the face mask.

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